

General System Requirements for Industrial Stand-Alone Software

Item	Recommended	Minimum
Processor	Intel Core 2 duo, AMD Athlon II or higher	Pentium IV
CPU Speed	2 Ghz or higher	1 Ghz
RAM	2 GB or higher	1 GB
Operating System	Windows® 7, Windows® 8, Windows® 10	Windows® 7, Windows® 8, Windows® 10
Hard Disk	5 GB free space on system disk	1 GB free space on system disk 2GB free space on data disk
CD	24X speed or faster	4X speed
Video Adapter	XGA [1024 x 768] 32 bit True Color Support for OpenGL	SVGA [1024 x 768] 16 bit High Color Support for OpenGL
Com Ports	1 Serial or USB port for each connected instrument	1 USB port
Internet / E-Mail	Recommended for data transfers, software updates, and technical support	Not required